Anytime Anywhere

Start Rental Business

- According to the lease agreement
 set the corresponding lease days and lighting times
- By setting the expired days of rental and the emitting shots to control whether the machine is running, more convenient for your leasing business.

VISIT MORE




Remote control

Patient database

- Remotely view User treatment data through mobile phones, and obtain effective information through data analysis.
- The patient's treatment parameters can be saved by the machine with one click and uploaded to the mobile patient database, and the patient data can be directly retrieved during the next treatment


LEARN MORE



CUSTOMIZATION IN 16 LANGUAGES

Built-in multiple languages to make the operation smoother more languages are being added...














Now anytime, anywhere Only need a mobile phone to control the machine



Super long warranty period

Contains main components



Remote Service



Maintenances ervice



2YEARS WARRANTY

























SKIN DETECTION FUNCTION

Support skin detection function, real-time tracking treatment effec



New features, stay tuned...



200x lens

Scalp hair follicle amplification and fine detection

50x lens

Detect scalp skin condition

PLHANDLE

Multifunctional beauty equipment

Optional spot size:

8*40mm²&10*50mm²

Importd UK lamp long shots time:

300,000--500,000 times

4filters:

420-1200nm

530-1200nm


560-1200nm

640-1200nm



300,000-500,000 times

IPL handle long shots time



Imported

UK LAMP

UK imported lamp
Can be used for a long time
glow many times
Ensure stable working status

Learn mons

FILTERS

- Permanent Hair Removal(HR)
- Skin Rejuvenation(SR)
- Use Vascular Removal(VR)
- Skin Tightening
- Wrinkle Removal
- Breast Lifting
- Acne Removal
- Pigment Removal
- 4 Age Spot Removal etc.









STABILIZE







TECHNICAL SPECIFICATION

Model Name DepiMEDiX_Android Diode-

Machine Power 2000W

IPL Energy Density 1-50J

IPL Power Supply 2000W

IPL Spot Size 8*40mm² & 10*50mm²

Frequency 1-10Hz

IPL Energy Density 1-50J

Electrical Requirements 110/230VAC,15/20A max,50/60Hz

IPL Spot Size 8*40mm²

IPL Filter 420-1200nm, 530-1200nm 560-1200nm, 640-1200nm

Display

Net Weight

Dimensions

Stand-by Working

13.3"&15.6" android screen

60kgs

50*60*110cm


Continuously 24 hours

MULTIPLE LANGUAGES AVAILABLE

4K 13.3 & 15.6 INCH

LCD SCREEN

High sensitiviy ,huranizeddesign .easier operaton









XXXXXXXXXXXXXXXX



Push data remotely



Energy

Piease erter a value from 1-300



Frequency

Please erter a value from 1-10



Online demo

Leam more

48 HUURS Spotgoods ship

Stable current management system

TAIWAN MW power







Denimedikandroid Family Series Shipped to anywhere in the world every day





Follow FDA | En Español

SEARCH

Home

Food

Drugs

Medical Device

Radiation-Emitting Products

Vaccines, Blood & Biologic

Animal & Veterinan

osmetics

obacco Products

Establishment Registration & Device Listing

P FDA Home O Medical Devices O Databases

Formal certification document issued by the U.S. Food and Drug Administration (FDA) to a product or device

that meets its certification requirements.

Device Class:

Regulation Number:

878,4810

fedical Specialty:

CONTROL OF TAXABLE CONTROLS

Registered Establishment Name:

BEIJING LASERTELL MEDICAL CO., LTD.

Registered Establishment Number:

3023343409

Premarket Submission Number:

K220381

Owner/Operator:

Beiling LaserTell Medical Co., Ltd

Owner/Operator Number:

10085308

Establishment Operations:

Manufacturer