







Treatment theory

The machine uses DPL precision skin rejuvenation technology that is narrow-spectrum photo skin rejuvenation (DPL), it refers to a new optical skin rejuvenation technology, different from traditional technology (such as IPL, INTENSE PULSE LIGHT), The new technology can stimulate 100nm narrow-spectrum pulsed light which also contains the absorption peaks of melanin and oxygen and hemoglobin, allowing accurate concentration of effective therapeutic energy, thus quickly and efficiently solve the facial freckles and telangiectasia (red face) problems.




The machine uses DPL precision skin rejuvenation technology that is narrow-spectrum photo skin rejuvenation (DPL), it refers to a new optical skin rejuvenation technology, different from traditional technology (such as IPL, INTENSE PULSE LIGHT), The new technology can stimulate 100nm narrow-spectrum pulsed light which also contains the absorption peaks of melanin and oxygen and hemoglobin, allowing accurate concentration of effective therapeutic energy, thus quickly and efficiently solve the facial freckles and telangiectasia (red face) problems.



EpiMED_PL™ rear side appearance




Handpiece introduction





Advantages of handpiece



Application



User interface







Handpiece selection

Mode selection

SHR/IPL operation interface


Specification

EpiMED_PL™
3000W
650-950nm(S-DPL)500-600nm;550-650nm;585-950nm 415-515nm(optional)
1-26J/cm2
16x50mm2 (SHR) 16x30mm2 (IPL/E-Light) 10x30mm2 (IPL/E-Light) Optional
10Hz
15ms
Single + multi pulse
TEC+AIR+WATER
20 hours continues working
12.4' true color LCD touch screen
100-240VAC, 20A max., 50/60Hz
20kgs
30*50*40cm

Comparison **Traditional IPL** Advantage **Image** EpiMED PL ltem 690-950nm(HR) 650-1200nm (HR) 500-600nm More narrow wavelength width make treatment result Wavelength 560-1200nm(SR) 550-650nm faster and better 585-950nm Identified handles balance counter. The distributor or No **Auto detection** Yes owner can rent or control easilier. Identified features make operation easier and treatment safer. This speeds up the operation and reduces treatment 0.3-1Hz 1-10Hz Frequency time. Compared to traditional IPL treatment, it takes 1/8th of **Fixed Operation** Fixed and sliding the time to treat the back, legs and arms. Water temperature sensor, water flow sensor and water No Water circulation level sensor make sure the machine works safely Yes without worrying about handle burns and burn patients. **Energy calibration** No Yes **Energy calibration settings in factory mode.** Treatment head connector with water lock allows safe **Normal quality** Connection Yes and easy exchange even when only single treatment

head use is permitted.

Customer feedback




Before & After



About us

Our factory



About us

Our certificates







Certificate of Compliance


















About us **Our service**

LaserTell has a strong R&D team, including electronics engineers (designers, AVR programming, circuit layout), construction engineers, UI designer engineers, QC control engineers, etc. We provide OEM/ODM designs and programs.)

About us

After-sale service

2 years warranty available with normal proper operation from factory (2xFree + 1xCharge). Labor and materials are free, but there are no consumable parts.

Why 2+1 year warranty? (6 simple reasons)

- 1. Almost all medical and cosmetic lasers typically have a lifespan of only 5 years, but the first 2 years warranty and full maintenance are key points. LaserTell offers a 2-year standard warranty and a 1-year extended warranty. (2xFree, 1xCharge)
- 2. Lifetime service (regardless of 3 years or 5 years or more...)
- 3. LaserTell has a full complement of certified technicians and engineers as well as local service providers to assist with all your repair and maintenance needs, much cheaper and faster than other suppliers. provide excellent service.
- 4.24 hours + 7 days after-sales service and technical support.
- 5. Service manual and safety instructions are available.
- 6. Comprehensive training for 5-7 days at the factory and clinic, including how to operate, maintain and adjust various treatments.

Contact

Beijing LaserTell Medical Co., Ltd.

Tel: +86-10-56370026 Fax: +86-10-56370036 Sales Director: Allen Wang



+86-186 1021 7037 (Mob./WeChat)





+86-186 1021 7037



lasertell0001



lasertell.info



www.LaserTell.com; www.shrlasermachine.com www.lasertell.jp



info@lasertell.com



