

CORE COMPETITIVENESS
**INDUSTRY LEADING
TECHNOLOGY**

DESIGNED BY LASERTELL IN UK, ASSEMBLED IN CHINA!

WE SHARE WE CARE



ANDROID

TRIO WAVELENGTH DIODE LASER+IPL




Anytime Anywhere

Start Rental Business

- According to the lease agreement
set the corresponding lease days and lighting times
- By setting the expired days of rental and the emitting shots to control
whether the machine is running, more convenient for your leasing business.

VISIT MORE




Remote control

Patient database

- Remotely view User treatment data through mobile phones, and obtain effective information through data analysis.
- The patient's treatment parameters can be saved by the machine with one click and uploaded to the mobile patient database, and the patient data can be directly retrieved during the next treatment

LEARN MORE ➤



CUSTOMIZATION IN 16 LANGUAGES

Built-in multiple languages to make the operation smoother
more languages are being added...



English



French



Russia



Danish



Dutch



Greek

**Now
anytime, anywhere
Only need a mobile phone
to control the machine**



Super long
warranty period

Contains main
components



Remote Service



Maintenances ervice



2YEARS WARRANTY

LOOKING FOR A BUSINESS PARTNER
"YOUR BRAND, OUR EXPERTISE"
OEM&ODM



Design the style



Confirm the design



Develop software



Confirm mould




Produce inner structure



Finish the sample



48 hours dispatch



Prepare spare frames



Stable business partnership

NEW

SKIN DETECTION FUNCTION

Support skin detection function, real-time tracking treatment effect



New features, stay tuned...



01 200x lens

Scalp hair follicle amplification and fine detection

02 50x lens

Detect scalp skin condition

50 Million shots

We always adhere to the use of real USA
COHERENT lasers, enhance customer trust

500w/800w/1000w/1200w/1600w/1800w/2000w

USA COHERENT laser


[Learn more](#)



FAC Technology
Fast Axis Collimator

Collimator Lens Focus Beam

Collimate the light path, reduce the divergence angle of the beam,
and turn the divergent light path into a parallel light path.
Smoother the edges of the laser beam.



After addition of a collimator lens, beam focus



Without optical shaping, the beam diverges rapidly

COOLING SYSTEM

ICE POINT PAINLESS

-16°





Camera



$\phi 6$



40J

Facial Tip Hair Removal
Removable Design

Handle

Multifunctional Detachable Design

IPL HANDLE

Multifunctional beauty equipment

Optional spot size :

8*40mm²&10*50mm²

Imported UK lamp long shots time :

300,000--500,000 times

4filters:

420-1200nm

530-1200nm

560-1200nm

640-1200nm



300,000-500,000 times

IPL handle long shots time



Imported

UK LAMP

UK imported lamp
Can be used for a long time
glow many times
Ensure stable working status

Learn mons

FILTERS

- 01 Permanent Hair Removal(HR)
- 02 Skin Rejuvenation(SR)
- 03 Vascular Removal(VR)
- 04 Skin Tightening
- 05 Wrinkle Removal
- 06 Breast Lifting
- 07 Acne Removal
- 08 Pigment Removal
- 09 Age Spot Removal etc.



OPEN UP, CONNECT ON

PAIR INSTANTLY ONCE OPENED,TURN OFF ONCE CLOSED





STABILIZE



SMOOTH




HANDLE SCREEN




HAIR REMOVAL PRINCIPLE


MAKE YOUR BEAUTY BUSINESS MORE EASIER



The laser selectively actson
The melanin in the hair folli-
cle,destroying
The-hair-growing area



Natural hair loss, preventregen-
eration, and achievethe purpose
ofhair remova



Natural hair loss, preventregen-
eration, and achievethe purpose
ofhair remova

755NM 808NM 1064NM WAVE ADVANTAGE

Speed 808nm: golden standard wavelength for all hair types

Alex 755nm: specific effective for blonde hair, fine hair

Long pulse 1064nm:specific effective for dark, tanned skin best quality ipl shr opt laser/ nd yag /alexandrite

Laser hair removal machine price





Straight hair

Curly hair

Fine hair


ALL HAIR TYPE APPLY

**MULTIPLE LANGUAGES
AVAILABLE**

4K 13.3 & 15.6 INCH

LCD SCREEN

High sensitivity ,humanized design .easier operatcn



TECHNICAL SPECIFICATION


Model Name	TruMEDi07_Android Diode+IPL	Frequency	1-10Hz
Machine Power	4500W	Electrical Requirements	110/230VAC,15/20A max,50/60Hz
Laser Stack Power	500W/800W/1000W/1200W/1600W/1800W/2000W	Display	13.3"&15.6"android screen
IPL Power Supply	2000W	Net Weight	60kgs
IPL Spot Size	8*40mm ² & 10*50mm ²	Dimensions	50*60*110cm
Diode EnergyDensity	1-150J/cm ²	Stand-by Working	Continuously 24 hours
IPL Energy Density	1-50J		
Diode Spot Size	16x24mm ² /φ 6mm		
IPL Spot Size	8*40mm ²		
IPL Filter	420-1200nm, 530-1200nm 560-1200nm, 640-1200nm		

THREE WAVELENGTHS


755NM For light skin

808NM For yellow skin

1064NM For darker skin




AMERICAN ORIGINAL




American Original Import

By Original American Coherent lasertechnology, the original bars are adopted in the laser generator to ensure that the output is stable and sufficient. The high quality bars guarantee a super long life span of the laser.



TECHNOLOGY OF BENEFIT

Acting on the superficial layer of the skin, combined With 20HZ slfhr® flashes, efficient, stable, painless.



WHAT IS THE DIODE ?

proprietary hair removal innovation is clinically proven method of laser hair removal that is virtually painless, and it is the safe laser hair removal method for all skin types, including dark and tanned skin. Its unique gradual heating and rapid/fast pulse delivery method effectively damages the hair follicle and hinders regrowth, while preventing injury to the surrounding tissue.

UPGRADE TO 4 WAVELENGTHS

755NM

The 755nm alexandrite performs hair removal on Caucasian to olive skin tones.

808NM

The 808nm gold standard in Permanent Hair Removal and suitable on all pigmented hair and skin types.

1064NM

The Yag 1064 wavelength is characterized by lower melanin absorption, specialized for darker skin types, making it a focused solution for darker skin types.

TECHNOLOGY OF BENEFIT

HIGH power TEC




ITALIAN water pump



MW power supply



The machine parts are selected from **international first-line brands**
only to create high-quality machines.



Push data remotely



Energy

Please enter a value from 1-300



Frequency

Please enter a value from 1-10



Cooling

Please enter a value fom 1-10

[Online demo](#)

[Leam more](#)




48 HOURS
SPOTGOODS SHIP

Stable current management system

TAIWAN MW power

▶ Learn about power






Medical grade water flow precision

Italian water pump

Learn about water pump ◀



**15 years factory, 2 years warranty, 102 countries,
60,000+ customers**

Diode LASER—Make Your Beauty Business More Easier




755nm



808nm




1064nm




1-2 TIMES

Hair reduced by 75% growth



3-4 TIMES

Less hair; Younger and smoother skin



4 TIMES ABOVE

Hair stops to grow back





TruMEDiO7_Android DepiMEDiX_Android Family Series

Shipped to anywhere in the world every day

The background of the advertisement shows several stacked shipping containers. Each container is wrapped in dark material and secured with black straps. Labels are visible on the containers, providing details about the contents and shipping instructions. The labels include text such as "Ctn No. : 4-18", "Ctn No. : 3-18", "Ctn No. : 1-18", "Ctn No. : 5-18", and "Ctn No. : 6-18". Some labels also specify "By Air" and "Handpiece/手柄" or "Water Circuit/Water Circuit".

CHOOSE THE FASTEST shipping method

Safe and secure, your trustworthy partner-Beijing LaserTell



Follow FDA | En Español

SEARCH

Home Food Drugs Medical Devices Radiation-Emitting Products Vaccines, Blood & Biologics Animal & Veterinary Cosmetics Tobacco Products

Establishment Registration & Device Listing

FDA Home Medical Devices Databases



Formal certification document issued by the U.S. Food and Drug Administration (FDA) to a product or device that meets its certification requirements.

Classification Name:	POWERED LASER SURGICAL INSTRUMENT
Device Class:	2
Regulation Number:	873.4810
Medical Specialty:	General & Plastic Surgery
Registered Establishment Name:	BEIJING LASERTELL MEDICAL CO., LTD.
Registered Establishment Number:	3023343409
Premarket Submission Number:	K220381
Owner/Operator:	Beijing LaserTell Medical Co., Ltd.
Owner/Operator Number:	10085308
Establishment Operations:	Manufacturer